Red Dot™ Electrodes**Intended Use:**

The 3M™ Red Dot™ Foam Monitoring Electrode 2228 and the 3M™ Red Dot™ Foam Monitoring Electrode with 4mm Adapter 2228BA are intended to be used by healthcare professionals for ECG monitoring of adults. These electrodes are disposable, intended for single use, and have been tested for up to 3 days wear.

Product Description:

- The 3M™ Red Dot™ Foam Monitoring Electrodes 2228 and 2228BA include a snap connector, sticky gel, and a silver/silver chloride sensing element. These electrodes contain a pressure sensitive adhesive that is recommended for diaphoretic patients.

Warning:

- The 3M™ Red Dot™ monitoring electrodes 2228 and 2228BA with a metal snap have not been tested for use during Magnetic Resonance Imaging (MRI) procedures. The use of conductive patient-connected devices and patient lead wires/electrodes in MRI procedures may result in serious patient burns. If a physician orders this device to be used during MRI procedures, the user should not allow loops in patient lead wires or allow the lead wires to come in contact with the patient's skin.

Precautions:

- To prevent dry-out, electrodes should be used within 30 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- Replace self-adhesive electrodes if they no longer stick firmly to the skin.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical area to avoid unwanted RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.
- To avoid low adhesion, poor trace quality, or risk of infection due to cross contamination, do not reuse.

Instructions For Use:

- To open bulk packs, locate notch on back side of pouch. Grasp pouch on both sides of notch and gently pull open.

Open Other Side

Open Here

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin impedance reduction should be part of the skin preparation. Use 3M™ Red Dot™ Trace Prep as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.
- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- For electrode removal, pull slowly and gently from the edge of the electrode, while folding the electrode back on itself and supporting the skin underneath.

Shelf Life:

For shelf life, refer to the expiration date that is printed on each package.

Disposal:

Dispose of contents/container in accordance with the local/regional/national/international regulations.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

| Symbol Title | Symbol | Description and Reference |
|---|--------|--|
| Authorized Representative in European Community | | Indicates the authorized representative in the European Community. ISO 15223, 5.1.2 |
| Batch code | | Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5 |
| Catalogue number | | Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6 |
| Caution | | Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4 |
| CE Mark | | Indicates conformity to European Union Medical Device Regulation or Directive. |
| Date of Manufacture | | Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3 |
| Do not re-use | | Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2 |
| Do not use if package is damaged or open | | Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8 |
| Green Dot | | Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe. |
| Importer | | Indicates the entity importing the medical device into the EU. |
| Manufacturer | | Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1 |
| Medical Device | | Indicates the item is a medical device. |
| Natural rubber latex is not present | | Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B |
| Rx Only | | Indicates the U.S. Federal Law restricts this device to sale by or on the order of a physician. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1). |
| Use-by date | | Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4 |

For more information see, HCBGRegulatory.3M.com

Thank you for choosing 3M™ Red Dot™ Electrodes. If you are not completely satisfied or if you have questions or comments about this product, in the USA, please call the 3M Health Care Customer Helpline at 1-800-228-3957.

Outside of the USA, please contact your 3M Distributor or your 3M Sales Representative.

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3M and Red Dot are trademarks of 3M.
© 2019, 3M. All rights reserved.
GmbH**