## Intended Use:

The 3M™ Red Dot™ Monitoring Electrode 2200 series\* are intended to be used by healthcare professionals for ECG monitoring of adults. These electrodes are disposable, intended for single use, and have been tested for up to 3 days wear. The 2245-50 and 2248-50 electrodes can also be used for ECG monitoring of pediatrics.

\*Note: This information does not apply to the 3M™ Red Dot™ ECG Monitoring Electrode 2268, 2228, or 2244.

## Product Description:

The information provided in this insert applies to the 3M™ Red Dot™ Monitoring Electrode 2200 series with snap connector, pressure sensitive adhesive, and silver/silver chloride sensing element.

## Warning:

This device with metal snap connector has not been tested for use during Magnetic Resonance Imaging (MRI) procedures.

The use of conductive patient-connected devices and patient lead wires/electrodes in MRI procedures may result in serious patient burns. If a physician orders this device to be used during MRI procedures, the user should not allow loops in patient lead wires or allow the lead wires to come in contact with the patient's skin.

## Precautions:

- To prevent dry-out, electrodes should be used within 30 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry.
- Replace self-adhesive electrodes if they no longer stick firmly to the skin.
- To minimize skin irritation:
  - Avoid placing an electrode on an irritated skin site.
  - Do not abrade a skin site more than one time.
  - Do not prepare pediatric skin with alcohol or an abrader
  - Avoid removing electrodes frequently and/or reapplying to the same skin site.
  - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
  - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical area to avoid unwanted RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.
- To avoid low adhesion, poor trace quality, or risk of infection due to cross contamination, do not reuse.

## Instructions For Use:

- Select and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality (adult patients only), skin impedance reduction should be part of the skin preparation. Use 3M™ Red Dot™ Trace Prep or the electrode's built-in abrader as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.
- Remove the paper liner from the tab if applicable.
- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- To remove electrode, pull slowly and gently from the tab or edge of the electrode, while folding the electrode back on itself and supporting the skin underneath.

## Shelf Life:

For shelf life, refer to the expiration date that is printed on each package.

## Disposal:

Dispose of contents/container in accordance with the local/regional/national/international regulations.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

## Symbol Glossary

Symbol Title	Symbol	Description and Reference
Authorized Representative in European Community		Indicates the authorized representative in the European Community. ISO 15223, 5.1.2
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Do not use if package is damaged or open		Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Importer		Indicates the entity importing the medical device into the EU.
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Medical Device		Indicates the item is a medical device.
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Rx Only		Indicates the U.S. Federal Law restricts this device to sale by or on the order of a physician. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4

For more information see, [HCBGRegulatory.3M.com](http://HCBGRegulatory.3M.com)

Thank you for choosing 3M™ Red Dot™ electrodes. If you are not completely satisfied or if you have questions or comments about this product, in the USA, please call the 3M Health Care Customer Helpline at 1-800-228-3957. Outside of the USA, please contact your 3M Distributor or your 3M Sales Representative.

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3M and Red Dot are trademarks of 3M.
© 2019, 3M. All rights reserved.
**3M Deutschland GmbH**