Red Dot™ Electrodes

Intended Use:

The 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 series are intended to be used by healthcare professionals for ECG monitoring of adults. These electrodes are disposable, intended for single use, and have been tested for up to 5 days wear.

Product Description:

The 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 series include a snap connector, silver/silver chloride sensing element, and built-in abrader. They are radiolucent and contain a pressure sensitive adhesive that is recommended for diaphoretic patients.

Caution:

- The 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 Series were determined to be MR-Conditional via non-clinical testing. For information on the use of a 3M™ Red Dot™ monitoring electrode 2570 series during an MRI procedure please see below.

Precautions:

- To prevent dry-out, electrodes should be used within 45 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not use an electrode if the gel is dry or if the gel column is dislodged from the electrode.
- Replace self-adhesive electrodes if they no longer stick firmly to the skin.
- To minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid removing electrodes frequently and/or reapplying to the same skin site.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Assess electrode sites periodically.
- During surgical procedures, place ECG electrodes as far as possible from the electrosurgical area to avoid unwanted RF current flow through the electrode site. Otherwise burns could result at the ECG electrode site.
- To avoid low adhesion, poor trace quality, or risk of infection due to cross contamination, do not reuse.

Instructions For Use:

- To open bulk packs, locate notch on back side of pouch. Grasp pouch on both sides of notch and gently pull open.

Open Other Side

Open Here

- Odaberite i pripremite mjesto postavljanja EKSelect and prepare an electrode site according to your health care facility's protocol for ECG monitoring or diagnostic procedures.
- The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin.
- Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin impedance reduction should be part of the skin preparation. Use 3M™ Red Dot™ Trace Prep or the electrode's built-in abrader as part of the skin preparation.
- Attach the lead wire to the electrode.
- Remove electrode from the liner. Avoid touching the adhesive surface.
- Place the center of the electrode over the prepared skin site.
- Smooth the electrode down using a circular motion. Avoid pressing on the center of the electrode.
- For electrode removal, pull slowly and gently from the edge of the electrode, while folding the electrode back on itself and supporting the skin underneath.

Important Notes:

3M™ MR Conditional EKG/ECG Electrodes have not been tested for MRI-related issues with any MR Conditional monitoring system. Therefore, the acceptable use of 3M™ MR Conditional EKG/ECG Electrodes when connected to an MR Conditional monitoring system is unknown.

Disconnected 3M™ MR Conditional Electrodes can be worn safely during an MRI examination according to the conditions of use presented in the MRI Safety Information labeling. Any deviation from those specific conditions may result in injury (e.g., burn) to the patient.

MRI Safety Information

MR Conditional

The 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 series was determined to be MR Conditional. Non-clinical testing demonstrated that the 3M™ Red Dot™ Monitoring Electrode with Foam Tape and Sticky Gel 2570 series is MR Conditional. A patient with these electrodes can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the 3M EKG/ECG electrode is expected to produce a maximum temperature rise of 2.5°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the 3M EKG/ECG electrode extends approximately 3-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

Shelf Life:

For shelf life, refer to the expiration date that is printed on each package.

Dispose:

Dispose of contents/container in accordance with the local/regional/national/international regulations.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Authorized Representative in European Community		Indicates the authorized representative in the European Community. ISO 15223, 5.1.2
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223, 5.1.6
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Do not use if package is damaged or open		Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Importer		Indicates the entity importing the medical device into the EU.
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Medical Device		Indicates the item is a medical device.
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Rx Only		Indicates the U.S. Federal Law restricts this device to sale by or for use of a physician. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4

For more information see, HCBGRegulatory.3M.com

Thank you for questioning 3M™ Red Dot™ Electrodes. If you are not completely satisfied or if you have questions or comments about this product, in the USA, please call the 3M Health Care Customer Helpline at 1-800-228-3957. Outside of the USA, please contact your 3M Distributor or your 3M Sales Representative.

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