MED AIRE FLEX WAVE

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Table of contents

Preface	12
Declaration of conformity	12
Service life	12
Disposal	12
Identification plate	12
Reuse.....	12
Warranty	12
Introduction	13
Intended purpose	13
Indication contraindication.....	13
Scope of delivery.....	13
Safety instructions.....	13
EMC	14
Setting up the system.....	14
Explanation of functions	15
Startup.....	16
Hand check.....	16
Alternating pressure function.....	16
Static function	16
Autofirm	16
Seat.....	16
Air flow	17
CPR (Cardiopulmonary resuscitation)	17
Transport function	17
Power failure	17
Alarm function.....	17
Cleaning, disinfection and care	17
Reuse/patient change	18
Maintenance and service	18
Storage	18
Disposal.....	18
Troubleshooting	19
Technical data.....	19

Preface

Thank you for choosing a product from Drive DeVilbiss.

You will not be disappointed by the design, functionality and quality of this product.

Please read these operating instructions carefully prior to initial use. They contain important safety information and valuable tips for the proper use and care. If you have questions or need additional information, contact your Drive DeVilbiss specialist dealer who supplied your product.

These operating instructions contain all of the important information for adjustment and operation. The most current operating instructions are available as a PDF online or can be requested from customer service.

Enlarged formats for visually impaired persons can be displayed. Blind patients must be instructed by their attendant. Repairs as well as some adjustments require special technical training and must therefore be made by the medical supplies dealer.

Declaration of conformity

For the described product there is conformity according to the Regulation on medical devices (EU) 2017/745.

Service life

Our company assumes a service life of five years for this product as long as the product is used as intended and all maintenance and service requirements are observed.


This service life may be significantly exceeded if the product is carefully handled, repaired, maintained and used and there are no technical limits based on the further development of science and technology.

The service life may also be significantly shortened through extreme and improper use, however. Determination of the service life by our company does not represent any additional guarantee.

Disposal

If you no longer need the product, please contact your specialist dealer or bring the product to your local disposal center.

Identification plate

1 Article name | 2 Article number | 3 Importer | 4 Serial number | 5 Manufacturer | 6 EC-REP | 7 UDI code | 8 Max. load | 9 Follow the instructions for use | 10 Medical device | 11 Production date (year/month)

Reuse

This product can be reused.

Please note that all technical documents needed for safe handling must also be given to the new user.

The product must be checked by the specialist dealer, cleaned, and restored to optimum condition.

Always follow the manufacturer's instructions during preparation; these are made available on request.

Accessories are available from your specialist dealer or at www.drivedevilbiss.de

Warranty

The warranty period for the product supplied by us is 24 months from the date of purchase. Defective parts will be replaced free of charge if it can be shown that there is a defect in materials or manufacturing.

Parts subject to wear are generally not covered during the normal guarantee period unless the parts need repair or replacement which became necessary as a clear, direct consequence of a manufacturing or material error.

We reserve the right of error and modifications on the basis of technical improvements and changes in design.

EXCLUSION OF WARRANTY

Parts subject to wear include (among others): air filter, cover and fixation bands

Introduction

The Drive Medical alternating pressure system Med Aire FlexWave is specially developed for the prophylaxis and treatment of decubitus, up to and including grade 4 decubitus according to EPUAP.

This system replaces the existing mattress in your bed. Equipped with a cell-on-cell system and a high-performance unit with an automatic start-stop mechanism, this system is ideally suited for long-term therapy and facilitates home care.

The Med Aire FlexWave can be operated in alternating pressure mode as well as in static mode. Through an integrated seat function, the pressure in the cells can be adapted to a seated position.

A special autofirm function facilitates the mobilization of the patient on this system. The system is suitable for patients from 20 kg to a max. of 200 kg.

IMPORTANT! The Med Aire FlexWave does not replace regular repositioning of the patient!

Intended purpose

Aids against pressure sores serve to prevent and support (benefit) the treatment of pressure ulcers for individuals who are bedridden or sit permanently and are available on the market in the most varied designs and qualities. The MED AIRE FlexWave pump fills or vents the air cells of the MED AIRE FlexWave mattress alternately with ambient air, so that there is constant local pressure relief through free positioning and continuous pressure distribution. Even these products cannot replace the patient regularly changing position. The MED AIRE FlexWave is an automatic system. It has sensor technology, which enables continuous, automatic adjustment of the internal pressure of the aid when the patient changes position.

Indication | contraindication

Indication:

- Decubitus prophylaxis for patients at high to very high risk
- Pressure distribution and relief for immobile and partially mobile patients
- Decubitus therapy up to and including grade 4 according to EPUAP.

Static pressure-relief positioning is suitable for:

- Patients who only occasionally tolerate dynamic alternating pressure or who do not tolerate it at all
- Pain patients

Kontraindikation:

- in the case of a body weight under 20 kg and over 200 kg

- in the case of unstable fractures (particularly of the back, cervical and lumbar spine)

in the case of neurological diseases in which dynamic alternating pressure should not be used, only static pressure-relief positioning should be used.

To avoid uncertainties regarding the indication and contraindication, we recommend consulting the attending physician.

Scope of delivery

Please check the contents for damage and completeness. If one of the parts listed below is missing, please contact your specialist dealer.

1 x	Med Aire FlexWave unit
1 x	Med Aire FlexWave mattress replacement system (including cell-on-cell system, connection hose)
1 x	Breathable, water-repellent and quilted cover
1 x	Instructions for use

Safety instructions

- Keep the unit away from water and/or other liquids.
- If water or another liquid gets into the device, disconnect it from the mains immediately.
- Do not allow the system to come into contact with open flames or smoldering objects.
- Do not operate the unit in the vicinity of flammable substances or chemicals => Risk of explosion!!!
- The pump housing may only be opened by authorized technical personnel => Risk of electric shock!!!
- Avoid heavy strain and/or damage to the housing and/or power cord.
- If the device fuses are defective, please have the

device checked by an authorized specialist dealer.

- Use this product only for the purpose described in the instructions for use.
- Operate the device only with the specified voltage (230 V/50 Hz).
- The unit should not be covered => Risk of overheating!
- Connect the device only to the mattress intended for it.
- This device is safeguarded against electrical shock according to the specifications for type BF devices! Product protection class I
- Do not allow the system to come into contact with sharp objects (knives, scissors, cannulas, etc.).
- This device should not be disposed of with household waste. More information is available from local waste management companies.
- When using side rails on adjustable beds for disabled persons, a sufficient distance between the top edge of the positioning system and the upper end of the bed rail must be ensured. (220 mm according to DIN 1970:2000). If necessary, corresponding side rail risers should/must be attached to protect the patient.
- Every serious incident which occurs in connection with the product must be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is located.

EMC

This product meets the protection requirements listed in the council directive on the approximation of the laws of the member states on electromagnetic compatibility (EMC, IEC 60601-1-2).

However, mutual interference between electrical devices, particularly when mobile telephones are used, cannot be completely ruled out under certain conditions.

Setting up the system

The Med Aire FlexWave replaces the existing mattress in your standard or hospital bed.

It does not impair the functionality of the standard or hospital bed.

- Take the Med Aire FlexWave out of the packaging.
- Check the contents for completeness and any visible damage.

- Now remove the existing mattress from the standard or hospital bed if you have not already done so.
- Place the mattress replacement system in the existing standard or hospital bed. The foot end is marked with 2 feet.
- Secure the mattress replacement system using the available straps to the movable parts of the slatted frame to avoid impairing the functionality of the system and the adjustable slatted frame.
- Place the unit using the hooks provided at the foot end, on the floor, or on an even surface.

Now attach the connection hose to the unit **(photo 2)**.

Important! The connection hose must audibly click into place in order to prevent an unintentional separation between the unit and the mattress! Kinking or twisting the connection hose must be avoided.

Now the Med Aire FlexWave is ready for use.

Explanation of functions

Weight setting

User weight setting

Autofirm function

Max. increase of cell pressure in all cells

Static and alternating pressure

Adjustment of static pressure-relief positioning and dynamic alternating pressure.

Lock

Automatically enabled after 30 seconds of inactivity or by pressing the „Tastensperre/Lock“ button.

Alternating cycle

Adjustable: 10/15/20/25 minutes

Alarm "Off"

Acoustic alarm signal can be muted and remains muted until the button is pressed once again!

Seat function

Increase in the cell pressure in a seated position

Startup

Connect the plug to the unit and an appropriate power source (**photo 4**).

Switch on the unit. The on and off switch is located on the side of the device, above the power cord (**photo 5**).

The mattress is now in the filling phase while the LED „Low Pressure“ is lit (**photo 6**).


Adjust the weight of the patient using the (+) and (-) button on the weight adjustment scale*; if the „Low Pressure“ LED goes out and the „Normal Pressure“ LED is lit, the patient can be placed on the system.

***The weight adjustment scale indicated is a guideline and does not reflect the actual pressure in the air chambers. Always set to the nearest weight and check pressure by hand (see HANDCHECK section) and correct upwards or downwards, if necessary.**

Hand check

After each change in pressure and position, use your hand to check the proper pressure setting at the buttocks/hips. When correctly adjusted, a hand should be able to be easily inserted between the patient and a deflated cell (**photo 7**).

Alternating pressure function

All air chambers are evenly inflated and deflated in a 1:1 ratio. An alternating cycle involves in each case the inflating and deflating of an air chamber and takes approx. 10/15/20/25 minutes, depending on the setting. The first 3 air chambers in the head area remain static.

The „Alternating Pressure“ LED indicates that the alternating pressure function is enabled.

Static function

The Med Aire FlexWave has a static function as a standard feature.

This is switched on using the button  .

The „Static“ LED indicates that the static function is enabled.

The air chambers all maintain the same pressure and ensure adjustable pressure-relief positioning.

In static pressure-relief mode, the deepest point on the patient (buttock/hip area) should sink into the system to approximately halfway down the upper cell system (approx. 3-4 cm).

Note: Perform a hand check to verify; see hand check section.

Autofirm

This function facilitates the positioning and mobilization of the patient through a firm and stable surface. The pressure within the cells is increased to a maximum level. For reasons of safety, this function automatically switches off after 30 minutes!

Seat

In this case, the pressure in the air chambers is increased and adapted to the short-term seated position. The objective is to prevent the patient from „bottoming out.“

The seat function can be manually switched on and off. **This function may only be used for short periods of time.**

For safety reasons, this function automatically switches off after 90 minutes.

Important!: The system is not suitable for long periods of sitting in bed!

Air flow

The alternating pressure system has nine air flow cells. This group of air chambers is microperforated at the edge and in this way, it ensures that the alternating pressure system is actively supplied with air.

CPR (Cardiopulmonary resuscitation)

Before resuscitation can/may be performed on the system, the system must be fully deflated.

For rapid deflation, pull down on the red tab printed with „CPR“ at the head end. The tab is on the same side as the hose system. The valves open and air escapes from all of the cells.

To supply the system with air once again, close the CPR valve again by pressing it firmly and securely fasten the tab marked „CPR“ once again (**photo 8**).

Transport function

Remove the hose system from the unit and promptly close it with the attached cap. Now the mattress is in transport mode. To return to normal operation, remove the cap on the connection hose and connect it back to the unit and start the unit.

Note: When the Med Aire FlexWave is in transport mode, it is possible to leave patients on the system for up to 12 hours. However, at this time, we recommend regularly monitoring the patient and the mattress replacement system.

Power failure

In the event of a power failure, set the system to transport mode. See section on transport function.

Alarm function

Low pressure

The „Low Pressure“ LED is lit whenever the unit detects too little air in the cells. An acoustic signal sounds at the same time if it has not been switched off beforehand by pressing the „Mute“ button.

Once the set pressure has been reached in the cells again, the display goes out and the acoustic signal stops. The „Normal Pressure“ LED is lit.

Important: During the filling phase, the „Low Pressure“ LED is lit; once the pressure in the air chambers has reached the set user weight, the „Low Pressure“ LED goes out.

System failure/service display

If the unit detects an internal error, the „System Failure“ LED is continuously lit and an acoustic signal additionally sounds. This can be muted using the „Mute“ button. The „System Failure“ LED will still be lit. Please notify your specialist dealer.

Power failure

If there is a power failure or if the unit is disconnected from the power supply, the „Power Failure“ LED begins to blink and an acoustic signal sounds.

This alarm stays on until the unit is reconnected to the mains or the On/Off switch is operated.

Cleaning, disinfection and care

Mattress overlay system

The mattress overlay system should not be saturated or soaked*, however it can be cleaned by hand. To do this, use a moist (not wet) cloth and a mild standard cleaning and disinfectant solution.

* There is a risk of liquid getting into the cell chambers and hoses.

Unit

- Before cleaning, disconnect the unit from the power supply. During cleaning, please ensure that no fluids/cleaning solutions come into contact with the plug.
- The unit should not be saturated or soaked.
- The unit may be cleaned by hand. To do this, use a moist (not wet) cloth and a mild household cleaning and disinfectant solution.
- Do not use any solutions which contain phenol or abrasive cleaners. This can corrode the surface of the unit.
- Allow the unit to dry thoroughly.
- Perform a function test after the device is thoroughly dry.

Cover

- The cover can be cleaned by hand. To do this, use a moist (not wet) cloth and a mild cleaning solution.
- You can detach the cover using the zipper and wash it at a max. temperature of 60° using a household detergent.
- Do not dry the mattress cover in the dryer.

Note!

If the alternating pressure system is machine-washed and -disinfected, a validated and automated method according to the Robert Koch Institute (RKI) and/or the Association for Applied Hygiene (VAH) is necessary.

Reuse/patient change

- In the case of reuse/patient change, the alternating pressure system must be disinfected and cleaned according to a validated method recognized by the Robert Koch Institute.
- During processing, take note of the medical device law, the „Guidelines for the Processing of Medical Devices“ as well as the „Guidelines for Hospital Hygiene and Infection Prevention“ of the Robert Koch Institute.

Maintenance and service

Maintenance

This product is subject to the Medical Device Law (MPG) and the Medical Devices Operator Ordinance (MPBetreibV). Internationally, the respective valid national regulations apply.

The functionality and proper condition must be checked before each use. The following points must be regularly checked:

- Check power cord and plug for wear and damage.
- Check the housing of the unit for damage.
- Check the alternation between air supply and ventilation in each air chamber circuit. The alternating cycle can be found under “Alternating pressure function”.
- Check all hoses and connections for leaks, kinks and breakage.
- Check all LEDs which must light up in their respective function.

Please contact your specialist dealer if you identify any damage or defects. During all visual inspections and maintenance, the unit must be disconnected from the power supply. Testing and assessments must be performed and documented by suitable technical staff in accordance with DIN EN 62353.

Service

Only original replacement parts and accessories from Drive Medical may be used. If the system develops a defect or fault, please contact a trusted specialist dealer. Med Aire Wave can be checked and repaired by an authorized specialist dealer or via your specialist dealer from Drive Medical.

Drive Medical recommends performing a technical safety check at least every 2 years.

Storage

- Remove the hose system from the compressor.
- Fully deflate the mattress overlay.
- After you have fully deflated the overlay, roll up the overlay, starting at the foot end. You can use the fastening strap to secure the rolled-up overlay to prevent the alternating pressure system from unrolling.

Disposal

Please contact your local waste management company regarding professional disposal. Please be aware that the alternating pressure overlay and the mattress cover may be carrying microbes. This can lead to transmission of infection. Please dispose of the overlay and mattress cover such that there is no risk to you and third parties.

Troubleshooting

Symptom	Remedial action
Unit does not work	<ul style="list-style-type: none"> • Check whether the plug is correctly inserted. Switch on the unit again. • If the LED on the On/Off switch does not light up, check the voltage of the outlet or insert the plug into another outlet. • If the LED on the On/Off switch lights up but the unit does not work, please contact your specialist dealer.
“Low pressure“ alarm	<ul style="list-style-type: none"> • Check whether the connection hose on the unit is correctly inserted. • The CPR valve must be closed.
Mattress does not pump up	<ul style="list-style-type: none"> • Ensure that all cells are correctly connected to each other and that the connection hoses are not twisted or kinked.
Patient bottoms out	<ul style="list-style-type: none"> • This indicates excessively low pressure in the air chambers. Increase the user weight on the unit and perform a hand check to check whether the patient is sinking down. • Check air filter (located above the on/off switch on the side of the unit); please clean or replace it. • If the unit continues to generate too little pressure, the unit must be replaced. Contact your specialist dealer
No alternating pressure takes place	<ul style="list-style-type: none"> • Ensure that the air supply hoses are not kinked or twisted. • The static function must be switched off. When the static function is switched on, no alternating pressure takes place. • Ensure that the autofirm and seat functions are switched off.

Technical data

UNIT		MATTRESS REPLACEMENT SYSTEM		ENVIRONMENTAL PARAMETERS	
Power supply	230 V / 50 hz; protection class 1; IPX0; type BF, AP/APG NO	Dimensions	200 x 88 x 20 cm	Temperatures	
				Operation	10 - 35°C
				Storage	-15 – 50°C
				Shipment	-15 – 70°C
Power consumption	max. 25 Watt	cells	20 cells of which 9 are air flow cells, 3 are static head cells	Relative humidity	
				Operation	20 – 80% => non-condensing
				Storage	10 – 90% => non-condensing
Weight	3,9 kg	Cell height	20 cm		
Dimensions	36.9 x 13.5 x 25 cm	Cover	Nylon/PU, breathable, washable cells Nylon/PU;		
Pressure range	10 - 50 mbar				
Alternating cycle	10/15/20/25 min.				

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Drive Medical GmbH & Co. KG
Emergo Europe B.V.
Airflo (Xiamen) Medical Co., Ltd.